
pyns Documentation

Release 0.5

Elliott L. Barcikowski

April 03, 2014

CONTENTS

1	pyns Package	1
1.1	Overview	1
1.2	Required Python Modules	1
1.3	Installation	2
1.4	Neuroshare API to pyns Package	2
1.5	Examples	2
1.6	Revision History	4
2	pyns.NSFile class	5
3	pyns.nsfile module	6
4	pyns.nsententity Module	8
5	pyns.nsparser Module	15
6	pyns.nsexceptions module	22
7	Indices and tables	23
	Python Module Index	24
	Index	25

PYNS PACKAGE

The pyns package native Python Neuroshare implementation.

1.1 Overview

The classes and functions provided by this interface do not exactly replicate the functions provided by the Neuroshare API. Instead of the functional Neuroshare API generally written in c, the pyns package was created using an object oriented programming style with a focus on having a form natural to the Python language.

This project was put in place by Ripple LLC., after noticing the growing use of the Python language in the academic and scientific communities.

1.2 Required Python Modules

The pyns package makes use of a few Python modules outside the default modules that would be present in any Python installation. These modules are standard for any scientific or computational use of the Python language and are likely to be present on your system.

These modules are:

- numpy - A collection of array manipulation objects and functions.

When data is returned from data in the pyns package, numpy arrays are used to package the data for easy use with matplotlib.

- matplotlib - A collection of plotting and analysis objects and functions.

These are not used explicitly in the pyns package, but they are used in all the provided examples.

- psutil - A module containing process utilities for Windows, Linux, and OS X

This module is used to check the available system memory. Cached data for spike waveforms has the potential to arbitrarily large, and a check is made to ensure enough physical memory is present on a system. However, files large enough to cause problems are unlikely.

On a Windows system the Python distribution Python(x, y) makes installing these modules and other useful analysis packages easy. More information may be found at: <http://www.pythonxy.com>.

If using Linux platforms, the needed packages depend slightly on the Linux flavor, though they represent the same libraries. On Ubuntu install the following packages: python-psutil, python-numpy, python-matplotlib. They may be installed using apt-get with the following command:

```
$ sudo apt-get install python-psutil python-numpy python-matplotlib
```

If using a Red Hat variant, install the packages `python-psutil`, `python-matplotlib`, and `numpy`. These may be installed using `yum` with the command:

```
$ sudo yum install python-psutil numpy python-matplotlib
```

1.3 Installation

To install `pyns` system-wide on your machine either run the command line `setup.py`. This command should install successfully on any platform. If run in `bash`, this would look like:

```
$ python setup.py install
```

If using a Windows system, a `msi`, produced by Python distutils, is also available.

1.4 Neuroshare API to pyns Package

For those familiar the standard Neuroshare API, the following is a mapping of the traditional Neuroshare functions to their `pyns` equivalent. The arguments and return values are generally not exactly the same due to the object oriented nature of the `pyns` package. See the docstrings of the classes and functions for more information. Also, some examples of general use of `pyns` are provided towards the end of this section.

- `ns_OpenFile` - `pyns` equivalent: `pyns.NSFile`
- `ns_GetFileInfo` - `pyns` equivalent: `pyns.NSFile.get_file_info()`
- `ns_GetEntityInfo` - `pyns` equivalent: `pyns.nsententity.Entity.get_entity_info()`
- `ns_GetSegmentInfo` - `pyns` equivalent: `pyns.nsententity.SegmentEntity.get_seg_source_info()`
- `ns_GetSegmentSourceInfo` - `pyns` equivalent: `pyns.nsententity.SegmentEntity.get_seg_source_info()`
- `ns_GetEventInfo` - `pyns` equivalent: `pyns.nsententity.EventEntity.get_event_info()`
- `ns_GetAnalogInfo` - `pyns` equivalent: `pyns.nsententity.AnalogEntity.get_analog_info()`
- `ns_GetNeuralInfo` - `pyns` equivalent: `pyns.nsententity.NeuralEntity.get_neural_info()`
- `ns_GetEventData` - `pyns` equivalent: `pyns.nsententity.EventEntity.get_event_info()`
- `ns_GetSegmentData` - `pyns` equivalent: `pyns.nsententity.SegmentEntity.get_segment_data()`
- `ns_GetAnalogData` - `pyns` equivalent: `pyns.nsententity.AnalogEntity.get_analog_data()`
- `ns_GetEventData` - `pyns` equivalent: `pyns.nsententity.EventEntity.get_event_data()`
- `ns_GetTimeByIndex` - `pyns` equivalent: `pyns.nsententity.Entity.get_time_by_index()`
- `ns_GetIndexByTime` - `pyns` equivalent: `pyns.nsententity.Entity.get_index_by_time()`

1.5 Examples

A few examples are provided to ease new users to `pyns` and Python to this package.

1.5.1 Simple Entity Dump

This example will open a .nev file, find the associated .nsx files and print out the result of the `pyns.nsentivity.Entity.get_entity_info()` function for each entity that was found.

```
"""This quick example opens a .nev file named 'test_data.nev'
and prints the ns_EntityInfo struct for each entity.
"""
from pyns.nsfile import NSFile

nsfile = NSFile('test_data.nev')
for entity in nsfile.get_entities():
    print entity.get_info()
```

1.5.2 Spike Plotting

Make use of matplotlib to plot the first electrode channel was found to have recorded spike events. This example makes use of the pyns functions `pyns.nsfile.NSFile.get_entities()` and `pyns.nsentivity.SegmentEntity.get_segment_data()`.

```
"""Example of using matplotlib to plot Neuroshare data, retrieved
with pyns. Makes use of a fictious nev file 'test_data.nev'.
"""
from matplotlib import pyplot
import numpy
import pyns

nsfile = pyns.NSFile('test_data.nev')
# get segment entities
segment_entities = [e for e in nsfile.get_entities() if e.entity_type == 3]
segment_entity = None
# pick the first segment entity with a non-zero item_count
for entity in segment_entities:
    if entity.item_count > 0:
        segment_entity = entity
        break
t = numpy.arange(0, 52, dtype=float)/30
# overlay all the spikes for this entity
for item in range(0, segment_entity.item_count):
    (ts, data, unit_id) = segment_entity.get_segment_data(item)
    pyplot.plot(t, data)
pyplot.show()
```

1.5.3 Plotting Continous Channels

Make use of matplotlib to plot the first continous channel. This example makes use of the pyns functions `pyns.nsfile.NSFile.get_entities()`, `pyns.nsentivity.AnalogEntity.get_analog_info()`, and `pyns.nsentivity.AnalogEntity.get_analog_data()`.

```
"""Example of plotting Analog entities (continuous channels) using pyns
and matplotlib. Makes use of a fictious nev file 'test_data.ns2'.
"""
from matplotlib import pyplot
import numpy
import pyns
```

```
nsfile = pyns.NSFile('test_data.ns2')
# get segment entities
analog_entities = [e for e in nsfile.get_entities() if e.entity_type == 2]
analog_entity = None
# pick the first analog entity with a non-zero item_count
for entity in analog_entities:
    if entity.item_count > 0:
        analog_entity = entity
        break
if analog_entity == None:
    exit(0)
analog_info = analog_entity.get_analog_info()
# only get the first 50k data points (other wise we could be reading
# lots and lots of data in memory at once
data = analog_entity.get_analog_data(0, 50000)
t = numpy.arange(0, len(data), dtype=float)/analog_info.sample_rate
pyplot.plot(t, data)
pyplot.xlabel('[s]')
pyplot.ylabel('[uV]')
pyplot.show()
```

1.6 Revision History

1.6.1 Version 0.5

- Feature: Added support for float based continuous files (.nf3)
- Feature: Added support for stimulation markers returned as segment entities

1.6.2 Version 0.4

- First publicly released pyns version
- Provides a complete version of Neuroshare API

PYNS.NSFILE CLASS

class `pyns.NSFile` (*filename*, *proc_single=False*)

General entry point to the pyns implementation of the Neuroshare API. This class loads all the NEV files associated with the file specified. The files are read and all the found entities are stored. This class provides the port of the `ns_GetFileInfo` function from the Neuroshare API

get_entities (*entity_type=None*)

Return iterator to entity list. If *entity_type* is specified return only entities of the desired type. *entity_type* should be one of the members of `EntityType`

Parameter: *entity_type* – member from static class `nsentity.EntityType`,
default=None

Returns: Iterator to entity list containing wanted entities

get_entity (*entity_index*)

Return the entity specified by entity index

get_entity_count ()

Utility function to return the number of entities found in all the files.

get_file_data (*ext*)

Utility function to get the `FileData` instance with the specified *file_type*. *file_type* should be one of `nev`, `ns?`

Parameter: *ext* – file extension

Returns: Associated `FileData` instance or `None` if it is not found

get_file_info ()

equivalent function to the Neuroshare `ns_GetFileInfo` function. Returns: `FileInfo` namedtuple with `ns_FILEINFO` data

get_time ()

return the datetime class corresponding to the origin time found in `nev` files. This corresponds to the starting time found from the system clock of the DAQ computer.

has_file_type (*file_type*)

Checks to see if `NEURALEV`, `NEURALSG`, or `NEURALCD` type files were found

PYNS.NSFILE MODULE

Classes related to handling .nev and .nsx files.

The class below nsfile.NSFile is generally the entry point for the pyns codes.

class pyns.nsfile.**FileData** (*parser*)

Internal data to be used by the File class that is needed to find desired NEV and NSX data. members:

parser – File parser from the pyns.parser module, one of NevParser, Nsx21Parser, Nsx22Parser

time_span – Holds the time span of data for this file in seconds

extension

returns the file extension for this file as stored in the parser class

file_type

Return the file type, one of NEURALEV, NEURALSG, NEURALCD

name

return the full path to this data file. We don't need to store this directly, we may just pass this from the parser object.

class pyns.nsfile.**FileInfo**

FileInfo(file_type, entity_count, timestamp_resolution, time_span, app_name, time_year, time_month, time_day, time_hour, time_min, time_sec, time_millisec, comment)

app_name

Alias for field number 4

comment

Alias for field number 12

entity_count

Alias for field number 1

file_type

Alias for field number 0

time_day

Alias for field number 7

time_hour

Alias for field number 8

time_millisec

Alias for field number 11

time_min

Alias for field number 9

time_month

Alias for field number 6

time_sec

Alias for field number 10

time_span

Alias for field number 3

time_year

Alias for field number 5

timestamp_resolution

Alias for field number 2

class `pyns.nsfile.NSFile` (*filename*, *proc_single=False*)

General entry point to the pyns implementation of the Neuroshare API. This class loads all the NEV files associated with the file specified. The files are read and all the found entities are stored. This class provides the port of the `ns_GetFileInfo` function from the Neuroshare API

get_entities (*entity_type=None*)

Return iterator to entity list. If *entity_type* is specified return only entities of the desired type. *entity_type* should be one of the members of `EntityType`

Parameter: *entity_type* – member from static class `nsentity.EntityType`,

default=None

Returns: Iterator to entity list containing wanted entities

get_entity (*entity_index*)

Return the entity specified by entity index

get_entity_count ()

Utility function to return the number of entities found in all the files.

get_file_data (*ext*)

Utility function to get the `FileData` instance with the specified *file_type*. *file_type* should be one of `nev`, `ns`?

Parameter: *ext* – file extension

Returns: Associated `FileData` instance or `None` if it is not found

get_file_info ()

equivalent function to the Neuroshare `ns_GetFileInfo` function. Returns: `FileInfo` namedtuple with `ns_FILEINFO` data

get_time ()

return the datetime class corresponding to the origin time found in nev files. This corresponds to the starting time found from the system clock of the DAQ computer.

has_file_type (*file_type*)

Checks to see if NEURALEV, NEURALSG, or NEURALCD type files were found

PYNS.NSENTITY MODULE

Collection of classes containing data and functions to hold entity information. These classes facility the reading of nev and nsx (2.1 and 2.2) files to access info functions as well as getting data, waveforms, and timestamps.

All of the Neuroshare API info and data functions are included in the following classes, but these classes extend this and may enable users more flexibility than that API.

class `pyns.nsententity.AnalogEntity` (*parser, electrode_id, units, channel_index, scale, electrode_label=None*)

Holds data and functions needed for analog entities found in NS[0-9] files. This file stores some additional information found in NSx2.2 style files that include a variety of information about probes in their CC headers. Looking in constructor documentation for more information about these variables.

get_analog_data (*start_index=0, index_count=None, use_scale=True*)

equivalent of the `ns_GetAnalogData`. returns the analog data for this entity starting from the `start_index` bin and returning the next `index_count` bins. If `index_count == None`, returns to the end of the analog waveform.

Parameters: `start_index` - first bin of returned waveform (default: 0) `index_count` - number of bins in returned waveform (default: until end) `use_scale` - True: scale to physical units, False: use ADC count.

Default: True

Returns: analog waveform for this entity. Result will length `index_count` or `start_index` to the end file.

get_analog_info ()

equivalent to the `ns_GetAnalogInfo` function from the Neuroshare API

get_extended_header ()

read the extended header (CC) for this entity and return it.

get_index_by_time (*time, flag=0*)

Return the analog index to the segment best matched by time.

Parameters: `time` - time in seconds from the start of data taking `flag` - flag to describe how to match files:

-1 - Return the data entry occurring before

or inclusive of time.

0 - Return the data entry occurring at or closest to time.

1 - Return the data entry occurring after or inclusive of time. (default)

Returns: best matched index

get_time_by_index (*index*)

Retrieves time from entity index.

Parameters: index - index for wanted item

Returns: timestamp in seconds from start of data taking

label

Return entity label to be consistent with Neuroshare DLL. In the case of NSx2.1 these labels are included in the CC header. In NSx2.1 we create this label based on the period field.

sample_freq

Return the sample frequency in Hz for this analog channel. This number is the same for all analog channels of the same type and depends on the period and the timestamp_resolution. Generally timestamp_resolution will be the same for all entities but the period (or how much a waveform is decimated will change). All of these things are stored in the parser.

class pyns.nsentivity.**AnalogInfo**

AnalogInfo(sample_rate, min_val, max_val, units, resolution, location_x, location_y, location_z, location_user, high_freq_corner, high_freq_order, high_filter_type, low_freq_corner, low_freq_order, low_filter_type, probe_info)

high_filter_type

Alias for field number 11

high_freq_corner

Alias for field number 9

high_freq_order

Alias for field number 10

location_user

Alias for field number 8

location_x

Alias for field number 5

location_y

Alias for field number 6

location_z

Alias for field number 7

low_filter_type

Alias for field number 14

low_freq_corner

Alias for field number 12

low_freq_order

Alias for field number 13

max_val

Alias for field number 2

min_val

Alias for field number 1

probe_info

Alias for field number 15

resolution

Alias for field number 4

sample_rate

Alias for field number 0

units

Alias for field number 3

class `pyns.nsentivity.Entity`(*parser, electrode_id*)

Base class for Neuroshare Entities. This is an abstract class that holds data and functions common to all neuroshare entities. Actual, instances will of one of the below derived classes (SegmentEntity, AnalogEntity, EventEntity, or NeuralEntity)

add_packet_data(*timestamp, packet_index*)

Add packet data to list. The length of this list should be the same as item_count

get_entity_info()

return the entity info for this entity

get_index_by_time(*time, flag=0*)

Equivalent to the Neuroshare function `ns_GetIndexByTime`. Return the segment index to the segment best matched by time. Note: This is overridden by each derived Entity class

Parameters: *time* - time in seconds from the start of data taking *flag* - flag to describe how to match files:

-1 - Return the data entry occurring before

or inclusive of time.

0 - Return the data entry occurring at or closest to time.

1 - Return the data entry occurring after or inclusive of time. (default)

Returns: best matched index

get_time_by_index(*index*)

Equivalent to the Neuroshare function `ns_GetTimeByIndex`. Returns the time in seconds from the start of data taking for the time corresponding to index. Note: This is overridden by each derived Entity class

Parameters: *index* - index of wanted item

Returns: time in seconds since the start of data taking

label

Return default entity label. This label is used in Segment and Neural entities. This is consistent with both Matlab code as well as Neuroshare DLL.

resize_packet_data()

Reset packet data to the number of items found

class `pyns.nsentivity.EntityInfo`

`EntityInfo(label, entity_type, item_count)`

entity_type

Alias for field number 1

item_count

Alias for field number 2

label

Alias for field number 0

class `pyns.nsentivity.EntityType`

static class that will to interface with entity types and provide a couple of utility functions like getting the equivalent strings

```

classmethod get_entity_id (entity_string)
    Return enum value from entity string, one of “analog”, “segment”, “neural”, “event”, or unknown

classmethod get_entity_string (entity_type)
    return string corresponding to entity type

class pyns.nsentivity.EventEntity (parser, mode)
    Holds data and function for digital event entities found in NEV files.

    get_event_data (packet_index)
        equivalent of the ns_GetEventData from the Nueroshare API

    get_event_info ()
        equivalent of the Neuroshare function ns_GetEventInfo

        Returns: EventInfo instance

    get_index_by_time (time, flag=0)
        Return the segment index to the segment best matched by time.

        Parameters: time - time in seconds from the start of data taking flag - flag to describe how to match files:

            -1 - Return the data entry occurring before
                or inclusive of time.

            0 - Return the data entry occurring at or closest to time.

            1 - Return the data entry occurring after or inclusive of time. (default)

        Returns: best matched index

    get_time_by_index (index)
        Retrieves time range from entity indexes

        Parameters: index - index for wanted item

        Returns: timestamp in seconds from start of data taking

    label
        Return label for EventEntity

class pyns.nsentivity.EventInfo
    EventInfo(event_type, min_data_length, max_data_length, csv_desc)

    csv_desc
        Alias for field number 3

    event_type
        Alias for field number 0

    max_data_length
        Alias for field number 2

    min_data_length
        Alias for field number 1

class pyns.nsentivity.EventType
    Static class to interface with digital event types

class pyns.nsentivity.NeuralEntity (parser, electrode_id, unit_class, segment_entity)
    Entity class for Neural style entities. One of these instances will be created for each class and each segment
    entity. Storing this data will allow a user to browser through spike waveforms organized by the segment classes.

```

get_index_by_time (*time*, *flag=0*)

Return index to the segment best matched by time.

Makes use of the held segment entity in the neural entity class. This called is wrapped to segment_entity.get_index_by_time

Parameters: time - time in seconds from the start of data taking flag - flag to describe how to match files:

-1 - Return the data entry occurring before

or inclusive of time.

0 - Return the data entry occurring at or closest to time.

1 - Return the data entry occurring after or inclusive of time. (default)

Returns: best matched index

get_neural_info ()

equivalent of the ns_GetNeuralInfo from Neuroshare API

get_time_by_index (*index*)

Retrieves time from entity index.

Makes use of the held segment entity that in the neural entity class. This call is wrapped to segment_entity.get_time_by_index

Parameters: index - index for wanted item

Returns: timestamp in seconds from start of data taking

label

Return label for this Neural Entity by using the contained segment entity

class pyns.nsentivity.**NeuralInfo**

NeuralInfo(source_entity_id, source_unit_id, probe_info)

probe_info

Alias for field number 2

source_entity_id

Alias for field number 0

source_unit_id

Alias for field number 1

class pyns.nsentivity.**SegSourceInfo**

SegSourceInfo(min_val, max_val, resolution, subsample_shift, location_x, location_y, location_z, location_user, high_freq_corner, high_freq_order, high_filter_type, low_freq_corner, low_freq_order, low_filter_type, probe_info)

high_filter_type

Alias for field number 10

high_freq_corner

Alias for field number 8

high_freq_order

Alias for field number 9

location_user

Alias for field number 7

location_x
Alias for field number 4

location_y
Alias for field number 5

location_z
Alias for field number 6

low_filter_type
Alias for field number 13

low_freq_corner
Alias for field number 11

low_freq_order
Alias for field number 12

max_val
Alias for field number 1

min_val
Alias for field number 0

probe_info
Alias for field number 14

resolution
Alias for field number 2

subsample_shift
Alias for field number 3

class `pyns.nsentivity.SegmentEntity` (*parser, electrode_id*)

Holds data and information for Segment or spike data found in NEV files. The timestamps of spikes as well as the position of waveforms in a data file are stored in the `packet_data` list

get_extended_headers ()
searches through the NEV file and finds all the extended headers associated with this `electrode_id`. This method may be slow. If that is the case, I will at the location of each extended header to the `SegmentEntity` class

get_index_by_time (*time, flag=0*)
Equivalent to the Neuroshare function `ns_GetIndexByTime`. Return the segment index to the segment best matched by time.

Parameters: `time` - time in seconds from the start of data taking `flag` - flag to describe how to match files:

-1 - Return the data entry occurring before

or inclusive of time.

0 - Return the data entry occurring at or closest to time.

1 - Return the data entry occurring after or inclusive of time. (default)

Returns: best matched index

get_seg_source_info ()
Equivalent to the `ns_GetSegSourceInfo` from the Neuroshare API. Returns information found in the NEUEVFLT header, an extended header in .nev files.

Returns: `SegSourceInfo` instance

get_segment_data (*index*)

return the segment waveform corresponding to index for this entity.

This function is the Neuroshare equivalent to the ns_GetSegmentData

Parameters: index - index of wanted item

Returns:

tuple - (timestamp, waveform, unit_id)

timestamp - time of this item in seconds from start of data taking

waveform - waveform found in data packet. The units may be uV (with Neural Data) or V (with stim data)

unit_class - classification of item (0-255)

get_segment_info ()

return the segment info struct for this entity

get_time_by_index (*index*)

Equivalent to the Neuroshare function ns_GetTimeByIndex. Returns the time in seconds from the start of data taking for the time corresponding to index.

Parameters: index - index of wanted item

Returns: time in seconds since the start of data taking

class pyns.nsententity.**SegmentInfo**

SegmentInfo(source_count, min_sample_count, max_sample_count, sample_rate, units)

max_sample_count

Alias for field number 2

min_sample_count

Alias for field number 1

sample_rate

Alias for field number 3

source_count

Alias for field number 0

units

Alias for field number 4

PYNS.NSPARSER MODULE

pyns.nsparser - parser header and data for .nev, .nsx2.1 and .nsx2.2 files.

The pyns.nsparser module provides the functionality for all the unlying reading of data for the pyns packet. Any use of the pyns packet will make use of these classes behind the scenes. Advanced uses may find using these classes directly as an easy way to get data out of .nev and .nsx files easily.

The heavy lifting is done with the classes NevParser, Nsx21Parser, and Nsx22Parser. These maybe interfaced with through the factory ParserFactory.

The header formats and sizes are found for all headers and data packets. These are found in the top of the files as ???????_FORMAT and ???????_SIZE constants. namedtuples are provided to return the header and packet data.

class pyns.nsparser.CC

CC(header_type, electrode_id, electrode_label, phys_conn, conn_pin, min_dig_value, max_dig_value, min_analog_value, max_analog_value, units, high_freq_corner, high_freq_order, high_filter_type, low_freq_corner, low_freq_order, low_filter_type)

conn_pin

Alias for field number 4

electrode_id

Alias for field number 1

electrode_label

Alias for field number 2

header_type

Alias for field number 0

high_filter_type

Alias for field number 12

high_freq_corner

Alias for field number 10

high_freq_order

Alias for field number 11

low_filter_type

Alias for field number 15

low_freq_corner

Alias for field number 13

low_freq_order

Alias for field number 14

max_analog_value
Alias for field number 8

max_dig_value
Alias for field number 6

min_analog_value
Alias for field number 7

min_dig_value
Alias for field number 5

phys_conn
Alias for field number 3

units
Alias for field number 9

class `pyns.nsparser.DIGLABEL`
`DIGLABEL(header_type, label, mode)`

header_type
Alias for field number 0

label
Alias for field number 1

mode
Alias for field number 2

class `pyns.nsparser.NEUEVFLT`
`NEUEVFLT(header_type, packet_id, high_freq_corner, high_freq_order, high_filter_type, low_freq_corner, low_freq_order, low_filter_type)`

header_type
Alias for field number 0

high_filter_type
Alias for field number 4

high_freq_corner
Alias for field number 2

high_freq_order
Alias for field number 3

low_filter_type
Alias for field number 7

low_freq_corner
Alias for field number 5

low_freq_order
Alias for field number 6

packet_id
Alias for field number 1

class `pyns.nsparser.NEUEVLBL`
`NEUEVLBL(header_type, packet_id, label)`

header_type
Alias for field number 0

label
Alias for field number 2

packet_id
Alias for field number 1

class `pyns.nsparser.NEUEVWAV`
`NEUEVWAV(header_type, packet_id, phys_conn, conn_pin, dig_factor, energy_thres, high_thres, low_thres, number_sorted_units, bytes_per_waveform, stim_amp_dig_factor)`

bytes_per_waveform
Alias for field number 9

conn_pin
Alias for field number 3

dig_factor
Alias for field number 4

energy_thres
Alias for field number 5

header_type
Alias for field number 0

high_thres
Alias for field number 6

low_thres
Alias for field number 7

number_sorted_units
Alias for field number 8

packet_id
Alias for field number 1

phys_conn
Alias for field number 2

stim_amp_dig_factor
Alias for field number 10

class `pyns.nsparser.NEURALCD`
`NEURALCD(header_type, maj_revision, min_revision, bytes_headers, label, comment, period, timestamp_resolution, time_origin, channel_count)`

bytes_headers
Alias for field number 3

channel_count
Alias for field number 9

comment
Alias for field number 5

header_type
Alias for field number 0

label
Alias for field number 4

maj_revision
Alias for field number 1

min_revision
Alias for field number 2

period
Alias for field number 6

time_origin
Alias for field number 8

timestamp_resolution
Alias for field number 7

class `pyns.nsparser.NEURALEV`
`NEURALEV(header_type, file_rev_major, file_rev_minor, file_flags, bytes_headers, bytes_data_packet, timestamp_resolution, sample_resolution, time_origin, application, comment, n_ext_headers)`

application
Alias for field number 9

bytes_data_packet
Alias for field number 5

bytes_headers
Alias for field number 4

comment
Alias for field number 10

file_flags
Alias for field number 3

file_rev_major
Alias for field number 1

file_rev_minor
Alias for field number 2

header_type
Alias for field number 0

n_ext_headers
Alias for field number 11

sample_resolution
Alias for field number 7

time_origin
Alias for field number 8

timestamp_resolution
Alias for field number 6

class `pyns.nsparser.NEURALSG`
`NEURALSG(header_type, label, period, channel_count, channel_id)`

channel_count
Alias for field number 3

channel_id
Alias for field number 4

header_type
Alias for field number 0

label
Alias for field number 1

period
Alias for field number 2

class `pyns.nsparser.NEVEvent`
`NEVEvent(timestamp, packet_id, reason, reserved, digital_input, input1, input2, input3, input4, input5)`

digital_input
Alias for field number 4

input1
Alias for field number 5

input2
Alias for field number 6

input3
Alias for field number 7

input4
Alias for field number 8

input5
Alias for field number 9

packet_id
Alias for field number 1

reason
Alias for field number 2

reserved
Alias for field number 3

timestamp
Alias for field number 0

`pyns.nsparser.NEVSegment`
alias of `Segment`

class `pyns.nsparser.NevParser` (*fid*)
Interface to .nev files. Also for easy reading of nev files and functions to read known basic and extended headers, as well as both segment data packets and event data packets. Member data is held to simply reading of headers and data packets.

file_type
Static function to return the 8 byte header associated with this file.

get_basic_header ()
Read the NEURALEV header from a nev file. input: fid id of NEURAL returns: struct containing nev header or packet data or None on failure

get_data_packet (*packet_index=None*)
Return the desired data packet. Parameter:

packet_index – index of desired data packet, default=None. if packet_index: seek to absolute path in file and return packet. else: return packet from current point in file.

Returns: NEVEvent or NEVSegment (depending on type of packet) containing data, None on failure.

get_data_packets()

Generator to loop over all data packets. Makes use the get_data_packets function

get_extended_header(*header_index=None*)

Return extended header for nev file from current position in file. This should be modified to be position independent with possibly a generator

get_extended_headers()

Generator to loop over all extended headers. Makes use of the get_extended_header function

get_packet_headers()

Defines an iterator that will return only the timestamps, electrode_id, and unit or reason for spike and digital channels respectively. This function will read in chunks in hopes to more efficiently (in terms of time) open NEV files.

class pyns.nsparser.**Nsx21Parser**(*fid*)

Interface to Nsx2.1 files.

Uses the Python struct module to read all the binary data found in the Nsx21 style files. Holds as member variables a file object and a small amount of header data to allow for easy retrieval.

file_type

Return 8 byte header for this file.

get_analog_data(*channel, start_index, index_count*)

Return the analog waveform for Nsx2.1 files. Returns data starting at the start_index bin and the next index_count bins. If the end of the file is reached before index_count return a waveform with as many bins as are found. scale provides the conversion from ADC counts to physical numbers. Parameters:

channel - index of the wanted electrode data start_index - first bin of electrode data to return

index - how many bins of the waveform to return

get_basic_header()

Return basic header for Nsx2.1 file. Returns NEURALSg instance or None with failure

scale

Returns the scale for Nsx2.1 files. As far as I can tell this is always 1 for NSx2.1

time_span

Return time_span of data in this file. Calculated from number of data points, period, and the clock speed.

timestamp_resolution

Return timestamp_resolution. For Nsx2.1 this quantity is not stored. It will always be 3000.0.

class pyns.nsparser.**Nsx22DataPacket**

Nsx22DataPacket(header, timestamp, n_data_points, data_points)

data_points

Alias for field number 3

header

Alias for field number 0

n_data_points

Alias for field number 2

timestamp

Alias for field number 1

class pyns.nsparser.**Nsx22Parser**(*fid*)

Interface to Nsx2.2 files.

Uses the Python struct module to read all the binary data found in the Nsx21 style files. Holds as member variables a file object and a small amount of header data to allow for easy retrieval.

file_type

Return 8 byte header for this file. Always NEURALCD

get_analog_data (*channel_index, start_index, index_count*)

Return the analog waveform for Nsx2.2 files. Returns data starting at the start_index bin and the next index_count bins. If the end of the file is reached before index_count return a waveform with as many bins as are found. Parameters:

channel - index of the wanted electrode data
start_index - first bin of electrode data to return
index - how many bins of the waveform to return

Returns: Requested waveform as a numpy.array

get_analog_packet (*channel_index, start_index, index_count*)

Generator to read a file in large chunks but return only the wanted channel and data from a time slice. Parameters:

channel_index - if specified, only return data for the wanted. If unspecified, returns all the data for a single time-slice

packet_index - return the wanted packet. If unspecified, returns as a yield. Otherwise,

get_basic_header ()

return the basic NEURALCD file header using the NEURALCD struct defined above.

get_data_packet (*packet_index=None*)

Return one full data packet. This will contain an array of all the digitized data for one moment in time. If packet_index is not specified we will read from the current moment in the file. Otherwise, we skip to the desired packet

get_data_packet_format (*index*)

calculate the number format for the index'th data packet using the number of data points and the channel count

get_extended_header (*header_index=None*)

Get the desired extended (CC) header. If header_index == None, then the next header is read from the current position of the file. If the header_index is >= 0 skip to that position in the file and return that CC header

get_extended_headers ()

generator to loop through extended headers. Makes use of the get_extended_header function.

n_data_points

Return the number of data points from each pause section

time_span

Return time_span of data in this file. Calculated from number of data points, period, and the clock speed

pyns.nsparser.ParserFactory (*filename*)

ParserFactory provides the interface to the Parser classes listed below and handles opening of and checking the type of the files. Based on the string found in the first few bytes, it returns correct class

PYNS.NSEXCEPTIONS MODULE

Collection of Error classes and return values replicating the Neuroshare return values.

Throughout the pyns code, instead of making use of return types, Python exceptions are raised when function fail. The exceptions in pyns classes will always be of the type `NeuroshareError`.

As often as possible, exceptions are raised with an error number consistent with the Neuroshare API. These are defined in the static class `NSReturnTypes`

class `pyns.nsexceptions.NSReturnTypes`
Neuroshare standard return values.

exception `pyns.nsexceptions.NeuroshareError`
Exception class to raise Neuroshare errors.

INDICES AND TABLES

- *genindex*
- *modindex*
- *search*

PYTHON MODULE INDEX

p

- `pyns`, [1](#)
- `pyns.nsententity`, [8](#)
- `pyns.nsexceptions`, [22](#)
- `pyns.nsfile`, [6](#)
- `pyns.nsparser`, [15](#)

INDEX

A

add_packet_data() (pyns.nsentivity.Entity method), 10
AnalogEntity (class in pyns.nsentivity), 8
AnalogInfo (class in pyns.nsentivity), 9
app_name (pyns.nsfile.FileInfo attribute), 6
application (pyns.nsparser.NEURALEV attribute), 18

B

bytes_data_packet (pyns.nsparser.NEURALEV attribute), 18
bytes_headers (pyns.nsparser.NEURALCD attribute), 17
bytes_headers (pyns.nsparser.NEURALEV attribute), 18
bytes_per_waveform (pyns.nsparser.NEUEVWAV attribute), 17

C

CC (class in pyns.nsparser), 15
channel_count (pyns.nsparser.NEURALCD attribute), 17
channel_count (pyns.nsparser.NEURALSG attribute), 18
channel_id (pyns.nsparser.NEURALSG attribute), 18
comment (pyns.nsfile.FileInfo attribute), 6
comment (pyns.nsparser.NEURALCD attribute), 17
comment (pyns.nsparser.NEURALEV attribute), 18
conn_pin (pyns.nsparser.CC attribute), 15
conn_pin (pyns.nsparser.NEUEVWAV attribute), 17
csv_desc (pyns.nsentivity.EventInfo attribute), 11

D

data_points (pyns.nsparser.Nsx22DataPacket attribute), 20
dig_factor (pyns.nsparser.NEUEVWAV attribute), 17
digital_input (pyns.nsparser.NEVEEvent attribute), 19
DIGLABEL (class in pyns.nsparser), 16

E

electrode_id (pyns.nsparser.CC attribute), 15
electrode_label (pyns.nsparser.CC attribute), 15
energy_thres (pyns.nsparser.NEUEVWAV attribute), 17
Entity (class in pyns.nsentivity), 10
entity_count (pyns.nsfile.FileInfo attribute), 6
entity_type (pyns.nsentivity.EntityInfo attribute), 10

EntityInfo (class in pyns.nsentivity), 10
EntityType (class in pyns.nsentivity), 10
event_type (pyns.nsentivity.EventInfo attribute), 11
EventEntity (class in pyns.nsentivity), 11
EventInfo (class in pyns.nsentivity), 11
EventType (class in pyns.nsentivity), 11
extension (pyns.nsfile.FileData attribute), 6

F

file_flags (pyns.nsparser.NEURALEV attribute), 18
file_rev_major (pyns.nsparser.NEURALEV attribute), 18
file_rev_minor (pyns.nsparser.NEURALEV attribute), 18
file_type (pyns.nsfile.FileData attribute), 6
file_type (pyns.nsfile.FileInfo attribute), 6
file_type (pyns.nsparser.NevParser attribute), 19
file_type (pyns.nsparser.Nsx21Parser attribute), 20
file_type (pyns.nsparser.Nsx22Parser attribute), 21
FileData (class in pyns.nsfile), 6
FileInfo (class in pyns.nsfile), 6

G

get_analog_data() (pyns.nsentivity.AnalogEntity method), 8
get_analog_data() (pyns.nsparser.Nsx21Parser method), 20
get_analog_data() (pyns.nsparser.Nsx22Parser method), 21
get_analog_info() (pyns.nsentivity.AnalogEntity method), 8
get_analog_packet() (pyns.nsparser.Nsx22Parser method), 21
get_basic_header() (pyns.nsparser.NevParser method), 19
get_basic_header() (pyns.nsparser.Nsx21Parser method), 20
get_basic_header() (pyns.nsparser.Nsx22Parser method), 21
get_data_packet() (pyns.nsparser.NevParser method), 19
get_data_packet() (pyns.nsparser.Nsx22Parser method), 21
get_data_packet_format() (pyns.nsparser.Nsx22Parser method), 21

- [get_data_packets\(\)](#) (pyns.nsparser.NevParser method), 19
[get_entities\(\)](#) (pyns.NSFile method), 5
[get_entities\(\)](#) (pyns.nsfile.NSFile method), 7
[get_entity\(\)](#) (pyns.NSFile method), 5
[get_entity\(\)](#) (pyns.nsfile.NSFile method), 7
[get_entity_count\(\)](#) (pyns.NSFile method), 5
[get_entity_count\(\)](#) (pyns.nsfile.NSFile method), 7
[get_entity_id\(\)](#) (pyns.nsentivity.EntityType class method), 10
[get_entity_info\(\)](#) (pyns.nsentivity.Entity method), 10
[get_entity_string\(\)](#) (pyns.nsentivity.EntityType class method), 11
[get_event_data\(\)](#) (pyns.nsentivity.EventEntity method), 11
[get_event_info\(\)](#) (pyns.nsentivity.EventEntity method), 11
[get_extended_header\(\)](#) (pyns.nsentivity.AnalogEntity method), 8
[get_extended_header\(\)](#) (pyns.nsparser.NevParser method), 20
[get_extended_header\(\)](#) (pyns.nsparser.Nsx22Parser method), 21
[get_extended_headers\(\)](#) (pyns.nsentivity.SegmentEntity method), 13
[get_extended_headers\(\)](#) (pyns.nsparser.NevParser method), 20
[get_extended_headers\(\)](#) (pyns.nsparser.Nsx22Parser method), 21
[get_file_data\(\)](#) (pyns.NSFile method), 5
[get_file_data\(\)](#) (pyns.nsfile.NSFile method), 7
[get_file_info\(\)](#) (pyns.NSFile method), 5
[get_file_info\(\)](#) (pyns.nsfile.NSFile method), 7
[get_index_by_time\(\)](#) (pyns.nsentivity.AnalogEntity method), 8
[get_index_by_time\(\)](#) (pyns.nsentivity.Entity method), 10
[get_index_by_time\(\)](#) (pyns.nsentivity.EventEntity method), 11
[get_index_by_time\(\)](#) (pyns.nsentivity.NeuralEntity method), 11
[get_index_by_time\(\)](#) (pyns.nsentivity.SegmentEntity method), 13
[get_neural_info\(\)](#) (pyns.nsentivity.NeuralEntity method), 12
[get_packet_headers\(\)](#) (pyns.nsparser.NevParser method), 20
[get_seg_source_info\(\)](#) (pyns.nsentivity.SegmentEntity method), 13
[get_segment_data\(\)](#) (pyns.nsentivity.SegmentEntity method), 13
[get_segment_info\(\)](#) (pyns.nsentivity.SegmentEntity method), 14
[get_time\(\)](#) (pyns.NSFile method), 5
[get_time\(\)](#) (pyns.nsfile.NSFile method), 7
[get_time_by_index\(\)](#) (pyns.nsentivity.AnalogEntity method), 8
[get_time_by_index\(\)](#) (pyns.nsentivity.Entity method), 10
[get_time_by_index\(\)](#) (pyns.nsentivity.EventEntity method), 11
[get_time_by_index\(\)](#) (pyns.nsentivity.NeuralEntity method), 12
[get_time_by_index\(\)](#) (pyns.nsentivity.SegmentEntity method), 14
- ## H
- [has_file_type\(\)](#) (pyns.NSFile method), 5
[has_file_type\(\)](#) (pyns.nsfile.NSFile method), 7
[header](#) (pyns.nsparser.Nsx22DataPacket attribute), 20
[header_type](#) (pyns.nsparser.CC attribute), 15
[header_type](#) (pyns.nsparser.DIGLABEL attribute), 16
[header_type](#) (pyns.nsparser.NEUEVFLT attribute), 16
[header_type](#) (pyns.nsparser.NEUEVLBL attribute), 16
[header_type](#) (pyns.nsparser.NEUEVWAV attribute), 17
[header_type](#) (pyns.nsparser.NEURALCD attribute), 17
[header_type](#) (pyns.nsparser.NEURALEV attribute), 18
[header_type](#) (pyns.nsparser.NEURALSG attribute), 18
[high_filter_type](#) (pyns.nsentivity.AnalogInfo attribute), 9
[high_filter_type](#) (pyns.nsentivity.SegSourceInfo attribute), 12
[high_filter_type](#) (pyns.nsparser.CC attribute), 15
[high_filter_type](#) (pyns.nsparser.NEUEVFLT attribute), 16
[high_freq_corner](#) (pyns.nsentivity.AnalogInfo attribute), 9
[high_freq_corner](#) (pyns.nsentivity.SegSourceInfo attribute), 12
[high_freq_corner](#) (pyns.nsparser.CC attribute), 15
[high_freq_corner](#) (pyns.nsparser.NEUEVFLT attribute), 16
[high_freq_order](#) (pyns.nsentivity.AnalogInfo attribute), 9
[high_freq_order](#) (pyns.nsentivity.SegSourceInfo attribute), 12
[high_freq_order](#) (pyns.nsparser.CC attribute), 15
[high_freq_order](#) (pyns.nsparser.NEUEVFLT attribute), 16
[high_thres](#) (pyns.nsparser.NEUEVWAV attribute), 17
- ## I
- [input1](#) (pyns.nsparser.NEVEvent attribute), 19
[input2](#) (pyns.nsparser.NEVEvent attribute), 19
[input3](#) (pyns.nsparser.NEVEvent attribute), 19
[input4](#) (pyns.nsparser.NEVEvent attribute), 19
[input5](#) (pyns.nsparser.NEVEvent attribute), 19
[item_count](#) (pyns.nsentivity.EntityInfo attribute), 10
- ## L
- [label](#) (pyns.nsentivity.AnalogEntity attribute), 9
[label](#) (pyns.nsentivity.Entity attribute), 10
[label](#) (pyns.nsentivity.EntityInfo attribute), 10
[label](#) (pyns.nsentivity.EventEntity attribute), 11
[label](#) (pyns.nsentivity.NeuralEntity attribute), 12
[label](#) (pyns.nsparser.DIGLABEL attribute), 16
[label](#) (pyns.nsparser.NEUEVLBL attribute), 16

label (pyns.nsparser.NEURALCD attribute), 17
label (pyns.nsparser.NEURALSG attribute), 18
location_user (pyns.nsententity.AnalogInfo attribute), 9
location_user (pyns.nsententity.SegSourceInfo attribute), 12
location_x (pyns.nsententity.AnalogInfo attribute), 9
location_x (pyns.nsententity.SegSourceInfo attribute), 12
location_y (pyns.nsententity.AnalogInfo attribute), 9
location_y (pyns.nsententity.SegSourceInfo attribute), 13
location_z (pyns.nsententity.AnalogInfo attribute), 9
location_z (pyns.nsententity.SegSourceInfo attribute), 13
low_filter_type (pyns.nsententity.AnalogInfo attribute), 9
low_filter_type (pyns.nsententity.SegSourceInfo attribute), 13
low_filter_type (pyns.nsparser.CC attribute), 15
low_filter_type (pyns.nsparser.NEUEVFLT attribute), 16
low_freq_corner (pyns.nsententity.AnalogInfo attribute), 9
low_freq_corner (pyns.nsententity.SegSourceInfo attribute), 13
low_freq_corner (pyns.nsparser.CC attribute), 15
low_freq_corner (pyns.nsparser.NEUEVFLT attribute), 16
low_freq_order (pyns.nsententity.AnalogInfo attribute), 9
low_freq_order (pyns.nsententity.SegSourceInfo attribute), 13
low_freq_order (pyns.nsparser.CC attribute), 15
low_freq_order (pyns.nsparser.NEUEVFLT attribute), 16
low_thres (pyns.nsparser.NEUEVWAV attribute), 17

M

maj_revision (pyns.nsparser.NEURALCD attribute), 17
max_analog_value (pyns.nsparser.CC attribute), 15
max_data_length (pyns.nsententity.EventInfo attribute), 11
max_dig_value (pyns.nsparser.CC attribute), 16
max_sample_count (pyns.nsententity.SegmentInfo attribute), 14
max_val (pyns.nsententity.AnalogInfo attribute), 9
max_val (pyns.nsententity.SegSourceInfo attribute), 13
min_analog_value (pyns.nsparser.CC attribute), 16
min_data_length (pyns.nsententity.EventInfo attribute), 11
min_dig_value (pyns.nsparser.CC attribute), 16
min_revision (pyns.nsparser.NEURALCD attribute), 17
min_sample_count (pyns.nsententity.SegmentInfo attribute), 14
min_val (pyns.nsententity.AnalogInfo attribute), 9
min_val (pyns.nsententity.SegSourceInfo attribute), 13
mode (pyns.nsparser.DIGLABEL attribute), 16

N

n_data_points (pyns.nsparser.Nsx22DataPacket attribute), 20
n_data_points (pyns.nsparser.Nsx22Parser attribute), 21
n_ext_headers (pyns.nsparser.NEURALEV attribute), 18
name (pyns.nsfile.FileData attribute), 6
NEUEVFLT (class in pyns.nsparser), 16

NEUEVLBL (class in pyns.nsparser), 16
NEUEVWAV (class in pyns.nsparser), 17
NEURALCD (class in pyns.nsparser), 17
NeuralEntity (class in pyns.nsententity), 11
NEURALEV (class in pyns.nsparser), 18
NeuralInfo (class in pyns.nsententity), 12
NEURALSG (class in pyns.nsparser), 18
NeuroshareError, 22
NEVEEvent (class in pyns.nsparser), 19
NevParser (class in pyns.nsparser), 19
NEVSegment (in module pyns.nsparser), 19
NSFile (class in pyns), 5
NSFile (class in pyns.nsfile), 7
NSReturnTypes (class in pyns.nsexceptions), 22
Nsx21Parser (class in pyns.nsparser), 20
Nsx22DataPacket (class in pyns.nsparser), 20
Nsx22Parser (class in pyns.nsparser), 20
number_sorted_units (pyns.nsparser.NEUEVWAV attribute), 17

P

packet_id (pyns.nsparser.NEUEVFLT attribute), 16
packet_id (pyns.nsparser.NEUEVLBL attribute), 17
packet_id (pyns.nsparser.NEUEVWAV attribute), 17
packet_id (pyns.nsparser.NEVEEvent attribute), 19
ParserFactory() (in module pyns.nsparser), 21
period (pyns.nsparser.NEURALCD attribute), 18
period (pyns.nsparser.NEURALSG attribute), 19
phys_conn (pyns.nsparser.CC attribute), 16
phys_conn (pyns.nsparser.NEUEVWAV attribute), 17
probe_info (pyns.nsententity.AnalogInfo attribute), 9
probe_info (pyns.nsententity.NeuralInfo attribute), 12
probe_info (pyns.nsententity.SegSourceInfo attribute), 13
pyns (module), 1
pyns.nsententity (module), 8
pyns.nsexceptions (module), 22
pyns.nsfile (module), 6
pyns.nsparser (module), 15

R

reason (pyns.nsparser.NEVEEvent attribute), 19
reserved (pyns.nsparser.NEVEEvent attribute), 19
resize_packet_data() (pyns.nsententity.Entity method), 10
resolution (pyns.nsententity.AnalogInfo attribute), 9
resolution (pyns.nsententity.SegSourceInfo attribute), 13

S

sample_freq (pyns.nsententity.AnalogEntity attribute), 9
sample_rate (pyns.nsententity.AnalogInfo attribute), 9
sample_rate (pyns.nsententity.SegmentInfo attribute), 14
sample_resolution (pyns.nsparser.NEURALEV attribute), 18
scale (pyns.nsparser.Nsx21Parser attribute), 20
SegmentEntity (class in pyns.nsententity), 13

SegmentInfo (class in pyns.nsententity), 14
SegSourceInfo (class in pyns.nsententity), 12
source_count (pyns.nsententity.SegmentInfo attribute), 14
source_entity_id (pyns.nsententity.NeuralInfo attribute), 12
source_unit_id (pyns.nsententity.NeuralInfo attribute), 12
stim_amp_dig_factor (pyns.nsparser.NEUEVWAV attribute), 17
subsample_shift (pyns.nsententity.SegSourceInfo attribute), 13

T

time_day (pyns.nsfile.FileInfo attribute), 6
time_hour (pyns.nsfile.FileInfo attribute), 6
time_millisec (pyns.nsfile.FileInfo attribute), 6
time_min (pyns.nsfile.FileInfo attribute), 6
time_month (pyns.nsfile.FileInfo attribute), 6
time_origin (pyns.nsparser.NEURALCD attribute), 18
time_origin (pyns.nsparser.NEURALEV attribute), 18
time_sec (pyns.nsfile.FileInfo attribute), 7
time_span (pyns.nsfile.FileInfo attribute), 7
time_span (pyns.nsparser.Nsx21Parser attribute), 20
time_span (pyns.nsparser.Nsx22Parser attribute), 21
time_year (pyns.nsfile.FileInfo attribute), 7
timestamp (pyns.nsparser.NEVEvent attribute), 19
timestamp (pyns.nsparser.Nsx22DataPacket attribute), 20
timestamp_resolution (pyns.nsfile.FileInfo attribute), 7
timestamp_resolution (pyns.nsparser.NEURALCD attribute), 18
timestamp_resolution (pyns.nsparser.NEURALEV attribute), 18
timestamp_resolution (pyns.nsparser.Nsx21Parser attribute), 20

U

units (pyns.nsententity.AnalogInfo attribute), 10
units (pyns.nsententity.SegmentInfo attribute), 14
units (pyns.nsparser.CC attribute), 16